

FDA-Approved Rapid HIV Antibody Screening Test –Purchasing Details

February 28, 2011

	<u>Price Per Device*</u>	<u>External Controls</u>	<u># of Tests per Case</u>	<u>Catalog Numbers</u>	<u>Storage Temperature</u>	<u>Operating Temperature</u>	<u>Shelf Life of Test**</u>	<u>Shelf Life of Control**</u>	<u>Total Time Required to Conduct Test***</u>	<u>Window Period for Reading Results****</u>
OraQuick ADVANCE Rapid HIV-1/2 Antibody Test www.orasure.com	\$17.50	Sold Separately (\$25)	25 or 100	#1001-0079 (25 tests) #1001-0078 (100 tests) #1001-0077 (controls)	2-27°C (tests) 2-8°C (controls)	15-37°C	24 months	12 months 21 days (after opening)	<5 minutes (<10 min. for plasma) +20 min wait time	20-40 min
Uni-Gold Recombigen HIV www.unigoldhiv.com	\$15.75 \$8.00*	Sold Separately (\$26.25)	20	#1206506 (tests) #1206530 (controls)	2-27 °C (tests) 2-8 °C (controls)	15-27°C	12 months	12 months 1 month (after opening)	<5 minutes + 10 min wait time	10-12 min
Reveal G-3 Rapid HIV-1 Antibody Test www.reveal-hiv.com	\$14.00	Included	20 or 60	B1057-6 (20 test kits) B1057-7 (60 test kits)	2-30°C (tests) 2-8°C (controls)	15-27°C	12 months	12 months 7 days (after reconstituting)	3-5 Minutes No add. wait time	Result must be read immediately

▪ “Public health” price for public health programs that are recipients of CDC funds for expanded HIV testing

*Actual price may vary by purchasing agreements with manufacturers

** From date of manufacture, unless otherwise noted

***First time listed is estimated time required to set up test. The second time is the required wait time before reading results. Times listed exclude time needed to draw/obtain sample

****As measured from last step of testing process

Note: Trade names are for identification purposes only and do not imply endorsement. This information was compiled from package inserts and direct calls to manufacturers.



Prepared by: Kali Stanger, Frances Margolin, and Jennifer Reiter at HRET; Bernard Branson, Jill Clark, and Margaret Lampe at CDC

FDA-Approved Rapid HIV Antibody Screening Test –Purchasing Details

February 28, 2011

	<u>Price Per Device*</u>	<u>External Controls</u>	<u># of Tests per Case</u>	<u>Catalog Numbers</u>	<u>Storage Temperature</u>	<u>Operating Temperature</u>	<u>Shelf Life of Test**</u>	<u>Shelf Life of Control**</u>	<u>Total Time Required to Conduct Test***</u>	<u>Window Period for Reading Results****</u>
MultiSpot HIV-1/HIV-2 Rapid Test www.biorad.com	\$25.00	Included	50	#72269 (test kits)	2-8 °C Or 20-30 °C (tests and controls)	20-30°C	12 months (@ 2-8 °C) 3 months (@ 20-30 °C)	12 months	10-15 minutes No add. wait time	Can be read immediately or anytime up to 24 hours
Clearview HIV 1/2 STAT-PAK www.alere.com	\$25.00 \$8.00	Sold Separately (\$50/set)	20 per kit/48 kits per carton	92110 (tests) 92112 (controls)	8-30°C or 46-86°F	18-30°C 64-86°F	24 months	24 months	<5 minutes + 15 min wait time	15-20 minutes
Clearview COMPLETE HIV 1/2 www.alere.com	\$26.00 \$9.00	Sold Separately (\$50/set)	25 per kit/48 kits per carton	92111 (tests) 92112 (controls)	8-30°C or 46-86°F	18-30°C 64-86°F	24 months	24 months	<5 minutes +15 min wait time	15-20 minutes

• “Public health” price for public health programs that are recipients of CDC funds for expanded HIV testing

* Actual price may vary by purchasing agreements with manufacturers

** From date of manufacture, unless otherwise noted

*** First time listed is estimated time required to set up test. The second time is the required wait time before reading results. Times listed exclude time needed to draw/obtain sample

**** As measured from last step of testing process

Note: Trade names are for identification purposes only and do not imply endorsement. This information was compiled from package inserts and direct calls to manufacturers.